

# DNA extraction protocol (QIAGEN DNeasy96)

Chang-Yu Chang

2021-02-08

## Aims

Extract DNA from pellets in 96-well plates.

## Materials

- PBS 36 mL
- proteinase K 2.5 mL
- buffer AL 20 mL
- ethanol 20 mL
- buffer AW1 50 mL
- buffer AW2 50 mL
- buffer AE 10 mL
- DNeasy96 plate 1
- S-block or DW96 1
- PCR plate 1

## Procedures

- ☐ Thaw pellet plates at room temperature.
- ☐ Dissolve pellets in 180 $\mu$ L of PBS.
- ☐ Transfer the dissolved pellets to DNeasy microtube plates.
- ☐ Add 25 $\mu$ L proteinase K and 200 $\mu$ L Buffer AL.
- ☐ Seal the plate with caps and mix by vortexing. Incubate at **56C** for **30 min**.
- ☐ Add 200 $\mu$ L pure ethanol.
- ☐ Seal the plate with caps and mix thoroughly by vortexing.
- ☐ Assemble the DNeasy 96 plates on top of DW96s. Mark the DNeasy 96 plates.
- ☐ Use mP1000 set at 625 $\mu$ L to transfer the lysate of each sample from DW96s to each well of DNeasy96 plates (maximum 900 $\mu$ L).
- ☐ Seal DNeasy96 with AirPore film. Use vacuum to collect the flow-through.
- ☐ Add 500 $\mu$ L of Buffer AW1 to each sample.
- ☐ Seal DNeasy96 with AirPore film. Use vacuum to collect the flow-through.
- ☐ Add 500 $\mu$ L of Buffer AW2 to each sample.
- ☐ Seal DNeasy96 with AirPore film. Use vacuum to collect the flow-through.
- ☐ Label two PCR plates. Place DNeasy96s on these PCR plates.
- ☐ Add 100 $\mu$ L Buffer AE to each sample, and seal the DNeasy 96 plates with new AirPore Tape Sheets (provided).
- ☐ Incubate for 1 min at room temperature (15–25°C). Centrifuge for **2 min** at **6000 rpm**.
- ☐ Cover the PCR plates with aluminum foil and save it in -20C freezer.